



GRINDING & POLISHING MACHINE

QPOL 30

The Qpol 30 is a transportable grinder and polisher designed for manual operation.

PRODUCT ADVANTAGES

- | Portable grinding device
- | Control unit for 2 micromotors
- | Increased torque due to reduction gear

PERFECT GRINDING - CONSUMABLES BY QATM

Grinding is more reliable, safe and efficient when processed with the right ingredients: knowledge, care and the right choice of consumables and instruments.



GRINDING & POLISHING MACHINE QPOL 30

TECHNICAL DATA

Speed	1000 - 40,000 rpm, continuously adjustable, incl. display
Wheel diameter	30 mm
Power connection	230 V / 50-60 Hz
W x H x D	130 x 97 x 254 mm
Weight	~ 2.3 kg

www.qatm.com/qpol-30

ORDER DATA